




Rapid - Acting Surface Disinfectant for Medical Devices Sensitive to Alcohol




Wear suitable protective gloves, goggles and clothing.



Spray the surfaces to be disinfected from a distance of 30 cm and ensure that the surface is completely covered.



Wait for the recommended action time.



Always complete the disinfection and cleaning process with a clean, dry cloth by wiping in an "S" shape from clean area to dirty, from peripheral to center and from top to bottom.

Actoanid® Spray is an alcohol-free, ready-to-use medical device disinfectant developed for use on all surfaces, including alcohol-sensitive surfaces. It is easy to use, fast effective and the smell is not disturbing.

Offers different usage preferences with spray and foam nozzle. Used safely on sensitive surfaces such as incubator and ultrasound probe.

Areas of Application:

Suitable for cleaning and disinfecting surfaces of alcohol-sensitive medical devices used in hospitals, doctor offices and clinics.

Method of Use:

Spray from 30 cm distance, until it covers the entire surface of the medical devices. Wait for the recommended exposure time.

Chemical Composition:

100 g contains: 0.23 g Poly (Hexamethylene Biguanide) Hydrochloride, (PHMB)
0.28 g Quaternary ammonium compounds, (QAC) auxiliary substances.

Microbiological Activity	Test Method	Concentration	Effect Time
Bactericidal (incl. MRSA)	EN 13727 EN 16615	Ready-to-use	1 min.
Bactericidal	EN 13697	Ready-to-use	5 min.
Yeasticidal (C.albicans)	EN 13624 EN 16615	Ready-to-use	1 min.
Yeasticidal (C.albicans)	EN 13697	Ready-to-use	5 min.
Limited Virucidal (Vaccinia-, Rota-, Corona viruses)	EN 14476	Ready-to-use	1 min.

▼
For warnings and precautions, see the label and user manual.