



» HYDROGEN PEROXIDE BASE SURFACE DISINFECTANTS «

# ACTO<sup>®</sup> OXY FOAM

**Disinfectant and Cleaning Foam with Sporocidal Effect for the Surfaces of Medical Devices**

- ✓ Contains hydrogen peroxide
- ✓ Ready to use and in foam form
- ✓ Suitable to use on surfaces sensitive to alcohol
- ✓ Aldehyde and phenol free
- ✓ Has wide material compliance
- ✓ Has bactericidal, fungicidal, virucidal and sporocidal



**Acto<sup>®</sup> Oxy Foam** is used for the cleaning and disinfection of the surfaces of medical devices. It is ready to use and in foam form. It is Hydrogen Peroxide based, but does not leave any residue. Thanks to its fast drying effect, it does not leave stains on the surfaces. It is bactericidal, fungicidal, virucidal and sporocidal.



German Company



**ACTO<sup>®</sup>** | Pharma  
Cosmetic  
Disinfection  
Wound Care  
Food Supplement  
Active Chemical Technology Organisation

# ACTO<sup>®</sup> OXY FOAM



## Disinfectant and Cleaning Foam with Sporocidal Effect for the Surfaces of Medical Devices

### Usage Area:

It is suitable for cleaning and disinfection of medical device surfaces (stethoscopes, cables, connectors, incubators, probes, etc.) that need to be disinfected frequently.

### Usage:


Spray from almost 30 cm away until it covers the whole surface of medical device. Make sure that the surface is completely wet. After waiting for the action time of the product, wipe the surface from remote to close and clean to dirty in one movement with a clean cloth and let it dry. Alternatively, spray the undiluted product on a cloth, wipe the surface and let it dry. Wait for the action time of it to take effect.

### Chemical Composition:


Content of 100g Solution:

2.6% Hydrogen peroxide, ethanol and auxiliary substances


### Instructions For Use:




Wear suitable protective gloves, goggles and clothing.



Spray the surfaces to be disinfected from a distance of 30 cm and ensure that the surface is completely covered.



Wait for the recommended action time.



Always complete the disinfection and cleaning process with a clean, dry cloth by wiping in an "S" shape from clean area to dirty, from peripheral to center and from top to bottom.

Packaging		
Unit	Pieces in Box	Article-No.
500 ml	18	05.5710.500
750 ml	15	05.5710.750
1 L	12	05.5710.1
4500+500 ml	3	05.5710.45
5 L	4	05.5710.5

Microbiological Activity	Test Method	Concentration	Effect Time
Bactericidal	EN 13727 EN 14561	Ready-to-use	1 min.
Fungicidal	EN 13624 EN 14562	Ready-to-use	1 min.
Sporicidal	EN 13704 EN 14347	Ready-to-use	30 min.
Virucidal (Adeno, Polio, MNV)	EN 14476	Ready-to-use	1 min.

