


## Ready-to-Use Disinfection Wipes for Alcohol-Sensitive Surfaces



**Actoanid® Wipes** is used for cleaning and disinfection of the surfaces of non-alcohol resistant medical devices which are not suitable for contact with liquid. It is polyhexanide and quaternary ammonium combination based. It is a product suitable for cleaning large areas and in the form of wipes containing nonwovens (viscose/polyester). It does not contain aldehyde and phenol derivatives.

### Usage:

Remove the blood and other body fluids on the surface with an absorbent material before the procedure. Wipe the surface of the medical device with **Actoanid® Wipes** to use 3-4 wipes per square meter. Make sure that the wipes reach the entire surface of the device and there is no dry place left. Wait for the duration of the effect. Drying should not be done after

use. The packaging should be tightly closed after each use. Used wipes should not be reused.

Microbiological Activity	Test Method	Concentration	Contact Time
Bactericidal	EN 16615	Ready-to-use	1 minute
Fungicidal (Aspergillus brasiliensis, Candida albicans)	EN 16615	Ready-to-use	1 minute
Virucidal (Vaccinia)	EN 14476	Ready-to-use	1 minute

### Chemical Composition:

Content of 100 g Impregnation Solution;  
0.25 g g Poly (Hexamethylene Biguanide) Hydrochloride, 0.25 g Quaternary Ammonium Compounds and excipients.

▼  
For  
warnings  
and  
precautions,  
see the  
label  
and user  
manual.